


# 3-Aminopyrrolidines



Racemic Compounds

Chiral Compounds

# 3-Aminopyrrolidines



















3-Aminopyrrolidines make up important components of bioactive compounds, and are found in pharmaceuticals and agricultural chemicals. The 3-aminopyrrolidine skeleton exerts a great influence on the pharmaceutical activity of 4-quinolonecarboxylate antibacterial agents. Moreover, 3-aminopyrrolidines have an asymmetric carbon atom and the relationship between their absolute configurations and the pharmaceutical activity has been studied thoroughly. With its established technologies, TCI has started commercial production of chiral 3-aminopyrrolidines, and is able to produce it in large quantities.





























Because 3-aminopyrrolidines have two amino groups in the molecule, they can generate di-substituted compounds and as well as two mono-substituted compounds by reaction with electrophiles. Therefore, if 3-aminopyrrolidines are to be introduced in a specific manner, it is necessary to protect one amino group. TCI has a wide assortment of mono-protected 3-aminopyrrolidines, where one amino group is protected by common protecting groups such as benzyl group or Boc group.

## References

- Tokyo Kasei Kogyo Co., Ltd., Jpn. Kokai Tokkyo Koho JP02218664, **1990**.  
 Tokyo Kasei Kogyo Co., Ltd., Jpn. Kokai Tokkyo Koho JP63041453, **1988**.  
 Tokyo Kasei Kogyo Co., Ltd., Jpn. Kokai Tokkyo Koho JP63041452, **1988**.  
 Tokyo Kasei Kogyo Co., Ltd., Eur. Pat. Appl. EP 218,249, **1987**.  
 M. A. Youngman, J. J. McNally, T. W. Lovenberg, A. B. Reitz, N. M. Willard, D. H. Nepomuceno, S. J. Wilson, J. J. Crooke, D. Rosenthal, A. H. Vaidya, S. L. Dax, *J. Med. Chem.* **2000**, *43*, 346.  
 J. P. Sanchez, R. D. Gogliotti, J. M. Domagala, S. J. Gracheck, M. D. Huband, J. A. Sesnie, M. A. Cohen, M. A. Shapiro, *J. Med. Chem.* **1995**, *38*, 4478.  
 J. P. Sanchez, A. J. Bridges, R. Bucsh, J. M. Domagala, R. D. Gogliotti, S. E. Hagen, C. L. Heifetz, E. T. Joannides, J. C. Sesnie, M. A. Shapiro, D. L. Szotek, *J. Med. Chem.* **1992**, *35*, 361.  
 T. Rosen, D. T. W. Chu, I. M. Lico, P. B. Fernandes, L. Shen, S. Borodkin, A. G. Pernet, *J. Med. Chem.* **1988**, *31*, 1586.

## Racemic Compounds

<b>Racemic Compounds</b>				
	<b>A1108</b> 5g 25g	<b>A1113</b> 5g	<b>A1112</b> 5g 25g	
				
	3-Acetamidopyrrolidine [79286-74-1]	3-(N-Acetyl-N-ethyl-amino)pyrrolidine [115445-29-9]	3-(N-Acetyl-N-methyl-amino)pyrrolidine [79286-87-6]	
<b>A1104</b> 5g 25g	<b>A1055</b> 5g 25g	<b>B1355</b> 25g	<b>B1354</b> 25g	<b>B1357</b> 10g
				
3-Aminopyrrolidine [79286-79-6]	3-Aminopyrrolidine Dihydrochloride [103831-11-4]	1-Benzyl-3-acetamidopyrrolidine [28506-01-6]	1-Benzyl-3-aminopyrrolidine [18471-40-4]	1-Benzyl-3-(tert-butoxycarbonyl-amino)pyrrolidine [99735-30-5]
<b>B1559</b> 5g	<b>B1560</b> 5g	<b>B1558</b> 5g 25g	<b>B1356</b> 25g	<b>B1358</b> 5g 25g
				
1-Benzyl-3-(dimethylamino)pyrrolidine [69478-77-9]	1-Benzyl-3-(ethylamino)pyrrolidine [115445-21-1]	1-Benzyl-3-(methylamino)pyrrolidine [96568-35-3]	1-Benzyl-3-(trifluoroacetamido)pyrrolidine [115445-23-3]	3-(tert-Butoxycarbonyl-amino)pyrrolidine [99724-19-3]
<b>M1096</b> 5g 25g	<b>D1964</b> 25g	<b>D1857</b> 5g	<b>D1859</b> 10g 25g	<b>D1965</b> 10g
				
3-(N-tert-Butoxycarbonyl-N-methylamino)pyrrolidine [172478-00-1]	N,N'-Dibenzyl-3-aminopyrrolidine [108963-20-8]	3-(Diethylamino)pyrrolidine [20984-81-0]	3-(Dimethylamino)pyrrolidine [69478-75-7]	N,N'-Dimethyl-3-aminopyrrolidine [64021-83-6]

<b>E0431</b> 5g  3-(Ethylamino)pyrrolidine [111390-22-8]	<b>M1106</b> 5g 25g  3-(Methylamino)pyrrolidine [83030-08-4]	<b>T1247</b> 10g  3-(Trifluoroacetamido)- pyrrolidine Hydrochloride [84424-06-6]			
<h2>Chiral Compounds</h2>			<b>A1169</b> 5g  (3R)-(+)-3-Acetamido- pyrrolidine [131900-62-4]	<b>A1170</b> 5g  (3S)-(-)-3-Acetamido- pyrrolidine [114636-31-6]	<b>A1167</b> 5g 25g  (3R)-(+)-3-Amino- pyrrolidine [116183-82-5]
<b>A1168</b> 5g 25g  (3S)-(-)-3-Amino- pyrrolidine [128345-57-3]	<b>A1053</b> 5g 25g  (3R)-(-)-3-Aminopyrrolidine Dihydrochloride [116183-81-4]	<b>A1054</b> 5g 25g  (3S)-(+)-3-Aminopyrrolidine Dihydrochloride [116183-83-6]	<b>A1173</b> 10g 25g  (3R)-(-)-1-Benzyl- 3-aminopyrrolidine [114715-39-8]	<b>A1174</b> 10g 25g  (3S)-(+)-1-Benzyl- 3-aminopyrrolidine [114715-38-7]	
<b>B1932</b> 5g  (3R)-(+)-1-Benzyl-3-(tert-butoxy- carbonylamino)pyrrolidine [131878-23-4]	<b>B1580</b> 5g  (3R)-(-)-1-Benzyl- 3-(ethylamino)pyrrolidine [376591-05-8]	<b>B1581</b> 5g  (3S)-(+)-1-Benzyl- 3-(ethylamino)pyrrolidine [169750-99-6]	<b>B1582</b> 5g 25g  (3R)-(-)-1-Benzyl- 3-(methylamino)pyrrolidine [144043-17-4]	<b>B1583</b> 5g 25g  (3S)-(+)-1-Benzyl- 3-(methylamino)pyrrolidine [169749-99-9]	
<b>B2920</b> 1g 5g  (3R)-(+)-1-(tert-Butoxycarbonyl)- 3-aminopyrrolidine [147081-49-0]	<b>B2921</b> 1g 5g  (3S)-(-)-1-(tert-Butoxycarbonyl)- 3-aminopyrrolidine [147081-44-5]	<b>A1171</b> 5g 25g  (3R)-(+)-3-(tert-Butoxycarbonyl)- amino)pyrrolidine [122536-77-0]	<b>A1172</b> 5g 25g  (3S)-(-)-3-(tert-Butoxycarbonyl)- amino)pyrrolidine [122536-76-9]	<b>D2149</b> 1g 5g  (3R)-(+)-3-(Dimethyl- amino)pyrrolidine [132958-72-6]	
<b>D2193</b> 1g 5g  (3S)-(-)-3-(Dimethyl- amino)pyrrolidine [132883-44-4]	<b>E0433</b> 1g  (3R)-(+)-3-(Ethylamino)- pyrrolidine [381670-30-0]	<b>E0434</b> 1g  (3S)-(-)-3-(Ethylamino)- pyrrolidine [381670-31-1]	<b>M1107</b> 1g 5g  (3R)-(+)-3-(Methylamino)- pyrrolidine [139015-33-1]	<b>M1108</b> 1g 5g  (3S)-(-)-3-(Methylamino)- pyrrolidine [139015-32-0]	
<b>T1369</b> 1g 25g  (3R)-(+)-3-(Trifluoroacetamido)- pyrrolidine Hydrochloride [141043-16-5]	<b>T1366</b> 1g 25g  (3S)-(-)-3-(Trifluoroacetamido)- pyrrolidine Hydrochloride [132883-43-3]				

**TCI** Ordering and  
Customer Service

**TCI AMERICA**

Tel : 800-423-8616 • 503-283-1681  
Fax : 888-520-1075 • 503-283-1987  
E-mail : sales@tciamerica.com  
Website : www.tciamerica.com

**East Coast Office**

Tel : 781-239-7515  
Fax : 781-239-7514

**TCI EUROPE N.V.**

Tel : 00 800 46 73 86 67 • +32 (0)3 735 07 00  
Fax : +32 (0)3 735 07 01  
E-mail : sales@tcieurope.eu  
Website : www.tcieurope.eu

**TCI Deutschland GmbH**

Tel : +49 (0) 6196 998678-0  
Fax : +49 (0) 6196 998678-1  
E-mail : sales@tcideutschland.de  
Website : www.tcieurope.eu/de/

**Tokyo Chemical Industry UK Ltd.**

Tel : +44 (0)1865 784560  
Fax : +44 (0)1865 784561  
E-mail : sales@tci-uk.co.uk  
Website : www.tci-uk.co.uk

**TOKYO CHEMICAL INDUSTRY CO., LTD.**

Tel : +81-3-5640-8878  
Fax : +81-3-5640-8902  
E-mail : globalbusiness@tokyokasei.co.jp  
Website : www.tci-asiapacific.com

**TCI Chemicals (India) Pvt. Ltd.**

Tel : 044-4261 2444  
Fax : 044-4261 1065  
E-mail : sales@tci-india.com  
Website : www.tci-india.com

**梯希爱(上海)化成工业发展有限公司**

Tel : 800-988-0390 • 021-6712-1386  
Fax : 021-6712-1385  
E-mail : sales@tcishanghai.com.cn  
Website : www.tcishanghai.com.cn